

This aims to be a very short introduction to how to use the pyxrr x-ray reflectivity simulation tool.

1. Important Files and Folders

| Base Directory: | |
|---------------------------------------|--|
| pyxrr/ | Folder containing the subroutines and modules for reflectivity calculation, the database and lots more |
| pyxrr/core.py | Python wrapper for xrr.c routines |
| pyxrr/functions.py | Helper functions for pyxrr.py |
| pyxrr/wrap4leastsq.py | Wrapper for scipy's leastsq fitting routine |
| pyxrr/xrr.c | C sources of reflectivity calculation routine $R(\theta)$ |
| pyxrr/xrr_sc.c | Same as above but for single core processing |
| pyxrr/materials.sqlite | Database containing dispersion corrections f_1+f_2 , atomic properties and densities |
| pyxrr/locale | Folder containing translations |
| | |
| contrib/ | Folder containing scripts and examples that use pyxrr and make it usable (the GUI must be placed here) |
| contrib/fit_interactive.py | Outer example script for interactive least squares fitting |
| contrib/fit_auto.py | Example for automated least squares fitting |
| contrib/multilayer_gammaopt.py | Outer example script for searching the optimum gamma parameter |
| contrib/pyxrr_GUI_v0.6.pyw | Graphical User Interface (GUI) for the simplified use of pyxrr |
| contrib/samples/ | Folder containing samples files defining the multilayer (see below) |
| contrib/samples/measurements/ | Folder containing the measured data files |
| contrib/results/ | Folder containing saved projects/results/simulations |
| contrib/tmp/ | Temporary folder containing results of last iteration |
| | |

| | |
|----------|---------------------|
| setup.py | Installation script |
|----------|---------------------|

2. Multilayer Sample Files

Every Multilayer consists of a Substrate, Groups of Layers and the Ambience. The Groups can be stacked to a given number of periods. An example for a samplefile is given below.

If a line ends with a comma (,) or a backslash (\) the next line is included in the processing.

EXAMPLE.param

Ambience: code=H2O, name=Water, rho=1.0

Group: name=protection, sigma=5, periods=1, grad_d=0

Layer: code=B4C, d=100, name=B4C_protection, rho=2.5

Group: name=ML, sigma=3, periods=100, grad_d=0.01

Layer: code=Mo, d=6, name=Mo_absorber, rho=9.4, sigma=1.5

Layer: code=B4C, d=3, name=B4C_barrier1, rho=2.3, sigma=1.3

Layer: code=Si, d=15, name=Si_spacer, rho=2.3, sigma=1.6

Layer: code=B4C, d=5, name=B4C_barrier2, rho=2.3, sigma=1.4

Substrate: code=SiO2, name=Glass, rho=2.34, sigma=1.1

Measurement: file=measurements/mosi_hzg4.njc, x_axis=twotheta, energy=8.905, fit_range=0.2->inf, resolution=0.01, background=-4., scale=1.0, offset=0, weighing=statistical, pol=0.5

Measurement: file=measurements/mosi_hasylib.fio, x_axis=theta, energy=10.0, fit_range=0.2->10., resolution=0.001, background=-5., scale=1.0, offset=0, weighing=z, pol=0

Measurement: energy=10.0

Explanation:

Italic grey attributes are **optional**. Their default values are listed below.

Red parameters are new in this version.

| Keyword | Definition | Default value |
|------------|--|---------------|
| code | Chemical formula of the material. Case sensitive. F.e.: La0.25Ca0.75MnO3 | - |
| name | Name to distinguish later. | code |
| periods | Number of Periods. How many repeats of the current group? | 1 |
| grad_d | Linear gradient in period thickness in % growth per layer. | 0 |
| sigma | RMS roughness of the upper interface in Å. | 0 |
| rho | density of the layer in g/cm ³ . | from database |
| d | Layer thickness in Å. | - |
| | | |
| file | Path to the file containing measured data. Can be .njc, .raw, .x00, .val, .asc, .fio or simple 2-columned data | - |
| x_axis | specify the independent (x-axis) quantity given in the measurement file. Possible values: <i>theta twotheta qz_A qz_nm</i> | |
| | | |
| | | |
| energy | Photon energy of the measurement in keV | 8 |
| twotheta | 1 if x-values are 2θ, 0 if they are θ | 0 |
| fit_range | θ_min->θ_max , Interval that will be evaluated for fitting simulated data (same unit as x_axis) | 0->inf |
| rebin | Rebinning is necessary if the measured data is unequally spaced in theta. Possible values: <i>average</i> (combine several data points to maximum step width) <i>interpolate</i> (interpolate data to minimum step width, resulting in more data points) | interpolate |
| resolution | FWHM of a gaussian shaped device resolution in θ/°. | 0 |
| background | Common logarithm of background Reflectivity (log10(R_noise)) | -10 |

| | | |
|-----------|---|--------------|
| scale | Scaling parameter for shifting Reflectivity values. (For Primary beam intensity uncertainties) Should be adjusted by hand in the beginning. | 1 |
| offset | θ offset in measured data in $^{\circ}$. | 0 |
| weighting | Weighting of the errorfunction. Possible values: z = Weights are given in 3rd column of measurement file statistical = $1/\sigma$, where σ is standard deviation of the Intensities (or proportional to it) | no weighting |
| pol | Polarization state of incident x-rays. (0 = perpendicular, 1 = parallel, 0.5 = unpolarized) This has negligible influence for $\theta < 10^{\circ}$ | 0 |

3.Theoretical Model (short)

pyxrr employs the Fresnel equations to calculate the reflection coefficients at each interface.

The interference of the reflected waves is calculated using the transfer matrix method (Abeles formalism).

Interface roughness is taken account for by implementing Nevot and Croce's formula for error-function shaped interfaces.

See for example:

http://en.wikipedia.org/wiki/Transfer-matrix_method_%28optics%29